



FOXsv


PTA Balloon Catheter

High Rated
Burst Pressure

Low tip entry profile

Excellent balloon
rewrap

Mastering the small vessel challenge



Designed for challenging small vessel PTA procedures

High Rated Burst Pressure – up to 22 bar

- Enables better treatment of highly calcified lesions

Low entry profile – 4F introducer compatible


- Eases insertion through tight or occluded lesions

Dual lumen shaft design – rapid inflation/deflation

- Enables shorter procedures and reduces patient discomfort

Excellent balloon rewrap – minimal deflated balloon profile

- Enables better recrossing of lesions/stents
- Eases insertion and withdrawal through introducer sheaths



**Endovascular
Commitment**

Abbott
A Promise for Life



Mastering the small vessel challenge

Recommended Guide Wire: 0.018"

Product Information

90 cm Usable Shaft Length Ref. No.	135 cm Usable Shaft Length Ref. No.	150 cm Usable Shaft Length Ref. No.	Balloon Diameter [mm]	Balloon Length [mm]	Sheath Compatibility [F]	Nominal Pressure (NP) [atm]	Rated Burst Pressure (RBP) [atm]
82158-54-01	82157-54-01	-	4.0	15	4	6	20
82168-54-01	82167-54-01	-	5.0	15	4	6	20
82144-54-01	-	83143-54-01	2.0	20	4	12	22
82373-54-01	-	83372-54-01	2.5	20	4	12	22
82150-54-01	-	83149-54-01	3.0	20	4	8	22
82160-54-01	82159-54-01	-	4.0	20	4	6	20
82170-54-01	82169-54-01	-	5.0	20	4	6	20
82180-54-01	82179-54-01	-	6.0	20	4	6	17
82146-54-01	-	83145-54-01	2.0	30	4	12	22
82375-54-01	-	83374-54-01	2.5	30	4	12	22
82152-54-01	-	83151-54-01	3.0	30	4	8	20
82162-54-01	82161-54-01	-	4.0	30	4	6	20
82172-54-01	82171-54-01	-	5.0	30	4	6	20
82148-54-01	-	83147-54-01	2.0	40	4	12	22
82377-54-01	-	83376-54-01	2.5	40	4	12	22
82154-54-01	-	83153-54-01	3.0	40	4	8	20
82164-54-01	82163-54-01	-	4.0	40	4	6	20
82174-54-01	82173-54-01	-	5.0	40	4	6	20
82184-54-01	82183-54-01	-	6.0	40	4	6	17
82198-54-01	-	83197-54-01	2.0	60	4	12	22
82379-54-01	-	83378-54-01	2.5	60	4	12	22
82156-54-01	-	83155-54-01	3.0	60	4	8	20
82166-54-01	82165-54-01	-	4.0	60	4	6	18
82176-54-01	82175-54-01	-	5.0	60	4	6	17
82186-54-01	82185-54-01	-	6.0	60	4	6	17
82946-54-01	-	83964-54-01	2.0	80	4	12	22
82949-54-01	-	83967-54-01	2.5	80	4	12	22
82952-54-01	-	83970-54-01	3.0	80	4	8	20
82955-54-01	82973-54-01	-	4.0	80	4	6	16
82958-54-01	82976-54-01	-	5.0	80	4	6	16
82961-54-01	82979-54-01	-	6.0	80	4	6	15
82947-54-01	-	83965-54-01	2.0	100	4	12	22
82950-54-01	-	83968-54-01	2.5	100	4	12	22
82953-54-01	-	83971-54-01	3.0	100	4	8	20
82956-54-01	82974-54-01	-	4.0	100	4	6	16
82959-54-01	82977-54-01	-	5.0	100	4	6	16
82962-54-01	82980-54-01	-	6.0	100	4	6	15
82948-54-01	-	83966-54-01	2.0	120	4	12	22
82951-54-01	-	83969-54-01	2.5	120	4	12	22
82954-54-01	-	83972-54-01	3.0	120	4	8	20
82957-54-01	82975-54-01	-	4.0	120	4	6	16
82960-54-01	82978-54-01	-	5.0	120	4	6	16
82963-54-01	82981-54-01	-	6.0	120	4	6	15

FOXsv is a trademark of the Abbott Group of Companies. Submarine Plus and Amphirion are trademarks of Invatec. Savvy is a trademark of Cordis. Symmetry is a trademark of Boston Scientific.

[Contact your local representative for more information](#)

Abbott Vascular Australia:

Level 8, 15 Talavera Road, North Ryde, NSW 2113. Tel: 1800 550 939, Tel: +61 2 8879 2837, Fax: +61 2 8879 2805

Abbott Vascular New Zealand:

Level 3, 6 Kingdom Street, New Market, Auckland. Tel: 0800 827 285, Tel: +64 9 524 6580, Fax: +64 9 524 6584

Photo(s) on file at Abbott Vascular. All drawings are artist's representations only and should not be considered as an engineering drawing or photograph. Products intended for use by or under the direction of a physician. Prior to use, it is important to read the package insert thoroughly for instructions for use, warnings and potential complications associated with use of this device. Information contained herein is for distribution for Australia & New Zealand ONLY.

For more information, visit our website at abbottvascular.com

© 2011 Abbott. All rights reserved. ANZ 060076 10/11

