

0.014"

0.018"

0.035"




SPARTACORE 14 STEELCORE 18 SUPRA CORE 35

Guide Wires


Extra support for
lower limb procedures




SPARTACORE 14



HI-TORQUE SteelCore 18 and 18 LT



HI-TORQUE Supra Core 35



**Endovascular
Commitment**

Abbott
A Promise for Life



Providing strong support and flexibility
for successful crossing

Product Information

SPARTACORE 14

Ref. No.	Floppy - Tip Length [cm]	Wire Length [cm]
1005201*	5	130
1005204*	10	130
1005202	5	190
1005205	10	190
1005203	5	300
1005206	10	300

Torque devices sold separately and in packs of ten (Ref. No. 22215) *DOC guide wire extension compatible (Ref. No.22260)

STEELCORE 18

Ref. No.	Tip Shape	Wire Length [cm]
1003281*	Straight	190
1003282	Straight	300

Torque devices sold separately and in packs of ten (Ref. No. 22215) *DOC guide wire extension compatible

STEELCORE 18 LT

Ref. No.	Tip Shape	Wire Length [cm]
1007709*	Straight	190
1007709-J*	J	190
1007710	Straight	300
1007710-J	J	300

Torque devices sold separately and in packs of ten (Ref. No. 22215) *DOC guide wire extension compatible

SUPRA CORE 35

Ref. No.	Tip Shape	Wire Length [cm]
1002703	Straight	145
1002703-01	Straight	190
1002703-02	Straight	300

Torque devices sold separately and in packs of ten (Ref. No. 1003279)

SPARTACORE 14, STEELCORE 18, STEELCORE 18 LT, SUPRA CORE 35 and **MICROGLIDE**
are trademarks of the Abbott Group of Companies.

[Contact your local representative for more information](#)

Abbott Vascular Australia:

Level 8, 15 Talavera Road, North Ryde, NSW 2113. Tel: 1800 550 939, Tel: +61 2 8879 2837, Fax: +61 2 8879 2805

Abbott Vascular New Zealand:

Level 3, 6 Kingdom Street, New Market, Auckland. Tel: 0800 827 285, Tel: +64 9 524 6580, Fax: +64 9 524 6584

Photo(s) on file at Abbott Vascular. All drawings are artist's representations only and should not be considered as an engineering drawing or photograph. Products intended for use by or under the direction of a physician. Prior to use, it is important to read the package insert thoroughly for instructions for use, warnings and potential complications associated with use of this device. Information contained herein is for distribution for Australia & New Zealand ONLY.

For more information, visit our website at abbottvascular.com

© 2011 Abbott. All rights reserved. ANZ 060075 10/11

